

Development of highly automated Next Generation Sequencing IVD tests for Solid Tumours

*Mengchu Wu, *Maurice Chan, *Arseny Smirnov, *David Hong, *Charlie Lee, *Pramila Ariyaratne, *Fabianus Mulawadi, *Sarawut Wongphayak, *Sameer Phalke, *Elaine Chua, *Atreyee Saha, *Tatiana Ivanova, *Wen Huang and *Gerd Michel

*Vela Research Pte Ltd., Singapore


Introduction

- With the increased use of NGS (Next Generation Sequencing) in clinical settings, Vela Diagnostics has developed four automated NGS based In Vitro diagnostics oncology disease panels (*Sentosa*[®] SQ Melanoma Panel, *Sentosa*[®] SQ CRC Panel, *Sentosa*[®] SQ NSCLC Panel and *Sentosa*[®] SQ Thyroid Cancer Panel) to aid medical decision-making.
- The panels are based on target amplification of disease related genes harbouring actionable/ prognostication related somatic hotspot mutations.

Design: Target gene list of the panels

| Targets for <i>Sentosa</i> [®] SQ Melanoma Panel | Targets for <i>Sentosa</i> [®] SQ CRC Panel | Targets for <i>Sentosa</i> [®] SQ NSCLC Panel | Targets for <i>Sentosa</i> [®] SQ Thyroid Panel |
|---|--|--|--|
| BRAF | NRAS | NRAS | NRAS |
| NRAS | CTNNB1 | CTNNB1 | CTNNB1 |
| CDKN2A | PIK3CA | PIK3CA | PIK3CA |
| MAP2K1 | FGFR3 | FGFR3 | FGFR3 |
| FGFR3 | KIT | KIT | EGFR |
| AKT3 | EGFR | EGFR | BRAF |
| KIT | BRAF | BRAF | RET |
| PIK3CA | RET | RET | PTEN |
| GNAQ | PTEN | PTEN | KRAS |
| GNA11 | KRAS | KRAS | TP53 |
| | TP53 | TP53 | |
| Total Number of Targeted Mutations | | | |
| 127 | 121 | 113 | 105 |

Methods and Workflow



Novel System Control


System Control (SC)

95% WT plasmid + 5% Mutant plasmid

- Negative control
- Extraction control
- Workflow (positive) control
- Spike-in control (detection of inhibitors in sample)
- Sequencing control (detect 5% of mutation)

| FFPE sample | SC-only sample (coverage > 1000) | Target Amplicons (median coverage ≥ 1000) | Interpretation |
|-------------|----------------------------------|---|---|
| + / - | + | + | Successful run |
| + | + | - | Suboptimal DNA input |
| - | + | - | PCR inhibitors in sample |
| + / - | - | - | General workflow failure |
| + / - | - | + | Possible workflow failure, results are not reliable |
| + / - | - | + | Possible contamination |

Automated and Integrated Workflow



Verification of *Sentosa*[®] SQ Panels Exemplified by *Sentosa*[®] SQ CRC Panel

| | | |
|---|--|--|
| LoD1 (min. DNA input for successful run) | 10 runs tested using Horizon Ref. Material carrying 10 known mutations, at a range of input (0.1 to 185ng) | PPA ¹ (95%CI) of 100.00 (95.91-100)% at 5 and 10ng |
| LoD2 (min. variant frequency VF for mutation detection) | 4 runs using Hor. Ref. Material carrying 39 mutations (Horizon 4-Pool gDNA Ref. Material), 24 with VF at ~5%. Tested at fixed input (LoD1) 9 runs using synthetic plasmids carrying 88 targeted mutations spiked into normal human genomic DNA at 5% VF | PPA (95%CI) of 100.00(64.6-100)% for VF ~5%. NPA (95%CI) of 100.00 (84.6-100)% for VF >7% |
| Specificity and reactivity | Bioinformatics analysis of 5 representative runs Unmapped reads ≤ 15%. Off-target reads ≤ 10%. No false-positive calls | Bioinformatics criteria met. |
| Stability test | Ambient temperature, freeze thaw cycles and accelerated stability study | 6 months |
| Interference | FFPE samples tested in 4 runs spiked with hemoglobin or triglycerides | Tolerance: At least 500ug (hemoglobin) or 250ug (triglycerides) |
| Repeatability | 18 runs using FFPE Ref. Materials using 3 machines, 3 operators, 2 lots | PPA (95%CI) 100 (99.7-100)% NPA (95%CI) 99.9 (99.8-99.9)% |

¹ PPA, positive percent agreement; NPA, negative percent agreement

Assay Specifications

| Assay Specifications | |
|--------------------------|---|
| Samples per Run | 1 system control + 7 patient samples |
| Sample Input | One 10 µm FFPE section |
| Mean Coverage / Sample | ≥ 1000x |
| TAT | 2.0 days |
| LoD (DNA input) | ≥ 5ng |
| LoD (mutation detection) | ≥ 5% variant frequency for mutation detection |
| Regulatory Status | |
| CE-IVD | |
| Australia-IVD | |

Clinical Validation of *Sentosa*[®] SQ Melanoma, CRC, NSCLC and Thyroid Cancer Panels

- Sentosa*[®] SQ Melanoma Panel**'s clinical sensitivity and specificity has been determined for the most frequent and clinically important target mutations in Melanoma samples – BRAF V600E
- The reference method is the CE-IVD *Sentosa* SA BRAF PCR Test
- Specimen: Melanoma FFPE (10 µm, ≥ 10% tumor content)
- 3 lots of *Sentosa* SQ Melanoma Panel Kit used

| <i>Sentosa</i> [®] SQ Melanoma Panel | | |
|---|--------------|--------------|
| Unique samples | BRAF V600E + | BRAF V600E - |
| 126 | 40 | 86 |
| BRAF V600E by qPCR | | |
| NGS Mutation calling at V600E | Positive | Negative |
| Positive | 76 | 1 |
| Negative | 0 | 95 |

BRAF status determined by *Sentosa* SA BRAF PCR Test

- Clinical sensitivity: **100%** (95% CI: 95.19 – 100%)
- Clinical specificity: **98.96%** (95% CI: 94.33 – 99.82%)
- In addition to BRAF V600E, 23 other mutations were called
- 13 from 51 variant detected were not confirmed by Sanger sequencing
- Mutation percentage varied between 1.56% – 6.65%

for 12 of these discordant calls

- Sentosa*[®] SQ CRC Panel**'s clinical sensitivity and specificity has been determined for the most frequent and clinically important target mutations in CRC samples – KRAS G12A/C/D/R/S/V/G13D
- The reference method is the CE-IVD *Sentosa* SA KRAS PCR Test
- Specimen: CRC FFPE (10 µm, ≥ 10% tumor content)
- 3 lots of *Sentosa* SQ CRC Panel Kit used

| <i>Sentosa</i> [®] SQ CRC Panel | | |
|--|-----------------|-----------------|
| Unique samples | KRAS (G12/13) + | KRAS (G12/13) - |
| 168 | 62 | 106 |
| KRAS G12A/C/D/R/S/V/G13D by qPCR | | |
| NGS Mutation calling at G12/13 | Positive | Negative |
| Positive | 78 | 0 |
| Negative | 0 | 101 |

KRAS status determined by *Sentosa* SA KRAS PCR Test

- Clinical sensitivity: **100%** (95% CI: 95.31 – 100%)
- Clinical specificity: **100%** (95% CI: 96.34 – 100%)
- In addition to KRAS G12A/C/D/R/S/V/G13D, 83 other mutations were called and tested using Sanger sequencing
- No discordance was identified

- Sentosa*[®] SQ NSCLC Panel**'s clinical sensitivity and specificity has been determined for the most clinically important and frequently mutated gene in NSCLC samples: EGFR
- The reference method is Sanger Sequencing
- Specimen: NSCLC FFPE (10 µm, ≥ 10% tumor content)
- 3 lots of *Sentosa* SQ NSCLC Panel Kit used

| <i>Sentosa</i> [®] SQ NSCLC Panel | | |
|--|----------------|----------------|
| Unique samples | EGFR targets + | EGFR targets - |
| 120 | 20 | 100 |
| EGFR targets by Sanger Sequencing | | |
| NGS Mutation calling at EGFR | Positive | Negative |
| Positive | 73 | 0 |
| Negative | 0 | 158 |

EGFR status determined by Sanger Sequencing

- Clinical sensitivity: **100%** (95% CI: 95 – 100%)
- Clinical specificity: **100%** (95% CI: 97.63 – 100%)
- The concordance (95% CI) for targeted mutations other than EGFR was 100% (95%CI: 80.64%-100%)
- No discordance was identified

- Sentosa*[®] SQ Thyroid Cancer Panel**'s clinical sensitivity and specificity has been determined for BRAF V600E
- The reference method was Sanger Sequencing
- Specimen: Thyroid cancer FFPE (10 µm, ≥ 10% tumor content)
- 3 lots of *Sentosa* SQ Thyroid Cancer Panel Kit used

| <i>Sentosa</i> [®] SQ Thyroid Cancer Panel | | |
|---|--------------|--------------|
| Unique samples | BRAF V600E + | BRAF V600E - |
| 62 | 11 | 51 |
| BRAF V600E by Sanger Sequencing | | |
| NGS Mutation calling at V600E | Positive | Negative |
| Positive | 74 | 0 |
| Negative | 0 | 80 |

BRAF status determined by Sanger Sequencing

- Clinical sensitivity: **100%** (95% CI: 95.07%-100%)
- Clinical specificity: **100%** (95% CI: 95.42%-100%)
- In addition to BRAF V600E, other detected mutations were called and tested using Sanger sequencing
- No discordance was identified

Conclusion

- All the panels (*viz.* *Sentosa*[®] SQ Melanoma Panel, *Sentosa*[®] SQ CRC Panel, *Sentosa*[®] SQ NSCLC Panel and *Sentosa*[®] SQ Thyroid Cancer Panel) have been verified and validated with ≥ 95% analytical and clinical sensitivity-specificity.
- In conclusion: robust, minimal hands-on, workflow-automated oncology NGS IVD panels show great promise to contribute significantly to personalized cancer diagnosis and treatment.