



Microbial Signatures and Genomic Profiling in Tumor Samples using Next Generation Sequencing.


Mei Qi Yee¹, Yik Lim Kok¹ Lina Nur Saleh¹, Tong Zhang², Pramila Ariyaratne¹, Yingnan Yu², Olivia Scully², Darwin Tay², Weng Kit Lye¹, Shi Hui Ang¹, Sirin Lee¹, Trang Nguyen¹, Claire Tang¹, Belinda Ling¹, Thaddaeus Ong¹, Harry Suhardi¹, Min Ko Ko Aye¹, Aung Win Kyaw¹, Jason Fischer³, Eugene Wee¹, Ludovic Lacroix⁴ & Charlie Lee^{1,2}
1 Vela Research, Vela Diagnostics Pte Ltd, Singapore, 2 Vela Genomics, Vela Diagnostics Pte Ltd, Singapore, 3 Vela Diagnostics Inc, USA, 4 Institut Gustave Roussy, France

Introduction

Precision medicine guides treatment options based on patients' genomic information. The advent of next-generation sequencing (NGS) has enabled the rapid acquisition of accurate genomic data. We have developed and validated a complete workflow (OncoKey SL Plus) comprising DNA and RNA extraction from the same formalin-fixed paraffin-embedded (FFPE) samples, library preparation, target enrichment, and high-throughput sequencing. Sequencing data is subsequently analyzed with an in-house bioinformatics pipeline to detect the presence of oncogenic bacteria and viruses, as well as an essential (OncoKey SL 60 Plus) or comprehensive (OncoKey SL 525 Plus) set of host's DNA mutations and RNA fusions associated with oncogenesis. Our workflows can be performed manually and have been optimised on the automated liquid handler Sentosa® SX 101 to prevent cross-contamination, with a 5 day turnaround time from 16 FFPE samples to lab report generation.

Materials and Methods

Complete Workflow Including Reporting and Interpretation for Oncology



Key features of OncoKey SL60 Plus & OncoKey SL525 Plus Panel

- ✓ Automation of 16 FFPE samples from sample preparation to results with 5 hours hands-on time and 5 day of turnaround time. (Including secondary analysis and interpretation with Sentosa SQ Reporter and Vela Analytics softwares, respectively.)
- ✓ Employs Unique Dual Indexes (UDI) and 10bp Unique Molecular Identifier (UMI) to enhance sensitivity and accuracy.
- ✓ Limit of detection: 40ng of DNA and RNA
- ✓ Detect SNP, Indels, CNV, Splice variants, Fusions, Microsatellite Instability (MSI) and Oncogenic pathogens
- ✓ Provide Tumor Mutational Burden value (OncoKey SL 525 Plus Panel)

- Standard reference materials used.

Customized FFPE reference standards	Horizon -HD-C331B
FFPE Tumor Fusion RNA reference materials	SeraSeq-0710-0496
Hepatitis B Virus (HBV) in PLC/PRF/5 cell line	ATCC-CRL-8024
Human papillomavirus (HPV)	Seracare-2400-0161
Merkel Cell Polyomavirus (MCPyV) in MKL-1 cell line	Sigma-9111801
HHV8 in BC-1 cell line	ATCC-CRL-2230
EBV in Raji cell line	ATCC-CCL-86
HCV Genotype 1a	Boca Biolistics

- 16 Horizon reference standards were used in the validation of the automated OncoKey SL assays. Sequencing was carried out using illumina MiSeq™/NextSeq™ systems.

- 16 SeraSeq® reference materials were used to assess the effectiveness of the assay in identifying splice variants and fusions. The success rate of detection was determined.

- Clinical samples, HBV, HPV, MCPyV, HHV8, EBV and HCV samples were used to assess the effectiveness of the assay to detect Oncoviruses.

Results

The tables below showed the validation results of the automated assays. Each run using 16 FFPE reference standards with LoD 40ng were performed. Data are representative of 10 repeated runs and showed a 100% pass in QC criteria for 8 DNA samples (left) and 8 RNA samples (right) in both assays.

OncoKey SL60 Plus		OncoKey SL525 Plus	
QC Performance	DNA Sample Barcode	QC Performance	DNA Sample Barcode
Mean Sample Coverage	689.89 - 851.01 (>200x)	Mean Sample Coverage	937.25 - 1269.22 (>200x)
Coverage Uniformity	99.32% - 99.46% (>90% targets with >50x coverage)	Coverage Uniformity	99.07% - 99.13% (>90% targets with >50x coverage)
Base Quality	89.66% - 91.41% (>80% bases with >30 quality score)	Base Quality	87.86% - 89.33% (>80% bases with >30 quality score)
Pass-Filter Reads	91.47% - 92.08% (>80% reads)	Pass-Filter Reads	100% (>80% reads)
Control Fragment Reads	87894 - 102446 (>1000 reads)	Control Fragment Reads	95858 - 102448 (>1000 reads)
Pass/Fail	All Pass	Pass/Fail	All Pass
Total Pass/Run	100%, 95% confidence interval	Total Pass/Run	100%, 95% confidence interval

Table 1 - OncoKey SL 60 Plus & OncoKey SL 525 Plus Quality Control Performance using Horizon DNA reference standards.

OncoKey SL60 Plus		OncoKey SL525 Plus	
QC Performance	RNA Sample Barcode	QC Performance	RNA Sample Barcode
Read Counts	616122 - 921262 (>10000 reads)	Read Counts	8155110 - 12182874 (>10000 reads)
Base Quality	87.43% - 90.50% (>80% bases with >30 quality score)	Base Quality	89.39% - 90.51% (>80% bases with >30 quality score)
Pass-Filter Reads	90.82% - 91.89% (>80% reads)	Pass-Filter Reads	100% (>80% reads)
Control Fragment Reads	9205 - 11780 (>1000 reads)	Control Fragment Reads	62190 - 94539 (>1000 reads)
Pass/Fail	All Pass	Pass/Fail	All Pass
Total Pass/Run	100%, 95% confidence interval	Total Pass/Run	100%, 95% confidence interval

Table 2 - OncoKey SL 60 Plus & OncoKey SL 525 Plus Quality Control Performance using SeraSeq® Fusion RNA reference materials.

Results demonstrated that the mutation detection sensitivity at 5% variant frequency for OncoKey SL60 Plus and SL525 Plus were 100% and 99%, respectively (left, using DNA reference standards). The mutation detection sensitivity at 5% variant frequency for OncoKey SL60 Plus and SL525 Plus were 96.4% and 97.9%, respectively (right, using RNA reference materials).

OncoKey SL60 Plus		OncoKey SL525 Plus	
Gene	Expected Mutation	Gene	Expected Mutation
BRCA1	V600E	BRCA1	V600E
EGFR	L858R	EGFR	L858R
KRAS	G12S	KRAS	G12S
NRAS	Q61H	NRAS	Q61H
PIK3CA	H1047R	PIK3CA	H1047R
MTT	Y981V	MTT	Y981V
CNA1	Q209L	CNA1	Q209L
Total Mutations (40-50%)	12	Total Mutations (40-50%)	12
Detected per Sample (excluding V600E* and L858R*+G12S*)	12	Detected per Sample (excluding V600E* and L858R*+G12S*)	12
Total Mutations (40-50%)	96	Total Mutations (40-50%)	96
Detected per Run	100%, 95% confidence interval	Detected per Run	99%, 95% confidence interval

Table 3 - The mutation detection efficiency of OncoKey SL 60 Plus & OncoKey SL 525 Plus using Horizon DNA reference standards. Read count with <1 are indicated as non detectable. * These samples with mutations have an expected allelic frequency of <5%, thus, they are not included in the acceptance criteria. Nonetheless, both OncoKey SL60 Plus and OncoKey SL525 Plus assay detected them.


OncoKey SL60 Plus		OncoKey SL525 Plus	
Gene	Expected Mutation	Gene	Expected Mutation
BRCA1	V600E	BRCA1	V600E
EGFR	L858R	EGFR	L858R
KRAS	G12S	KRAS	G12S
NRAS	Q61H	NRAS	Q61H
PIK3CA	H1047R	PIK3CA	H1047R
MTT	Y981V	MTT	Y981V
CNA1	Q209L	CNA1	Q209L
Total Mutations (40-50%)	12	Total Mutations (40-50%)	12
Detected per Sample (excluding V600E* and L858R*+G12S*)	14	Detected per Sample (excluding V600E* and L858R*+G12S*)	17
Total Mutations (40-50%)	144	Total Mutations (40-50%)	141
Detected per Run	96.4%, 95% confidence interval	Detected per Run	97.9%, 95% confidence interval

Table 4 - The mutation detection efficiency of OncoKey SL 60 Plus & OncoKey SL 525 Plus using SeraSeq® Fusion RNA reference materials. Read count with <1 are indicated as non detectable. # These RNA fusions are not targeted in the OncoKey SL 60 Plus assay. * Probe targeting this specific RNA fusion were not included in the particular run.

Results showed assay ability to detect oncogenic viruses- HBV, MCPyV, HPV51, HPV18 & HPV16. Previously the presence and the copy number of the extracted virus was determined using ddPCR.

Virus	Detected
HBV	✓
MCPyV	✓
HPV16,18,51	✓
HHV8	✓
EBV	✓
HCV	✓

32 FFPE clinical samples were analysed and HPV16 and EBV were identified in 2 head and neck samples. The presence of HPV and EBV was confirmed using immunohistochemistry (provided by Institut Gustave Roussy). OncoKey specifically identified different gene mutation PIK3CA p.R88Q and BRCA2 p.G602fs in the samples and provided detailed information on targeted therapy.



Conclusion.

Identifying virus and bacteria as etiological agents of specific cancers aids our understanding of the interplay between the microbiome and cancer, diagnosis and treatment. Robust automation coupled with UMI-based hybridisation capture workflow and NGS, we demonstrated that OncoKey SL60 and SL525 Plus assays can accurately identify over 60 clinically relevant DNA and RNA variants in FFPE samples, as well as the microbial signatures in various tumor types. Hence, we introduced a cost-effective strategy that enables high-throughput and rapid detection with accuracy and precision for comprehensive microbial and genomic profiling.