A Multiplex Molecular Diagnostic Platform and Test for Identification of Staphylococcus species and Detection of the mecA gene Directly from Positive Blood Cultures

Abstract


Background: Staphylococcal blood stream infections are increasingly prevalent and cause high morbidity and mortality. An automated molecular diagnostic platform, Portrait Dx, consists of a disposable cartridge run by a small instrument. The system combines multiplex amplification with multiplex chip-based detection. Here we describe development of the Staph ID/R assay for simultaneous identification of 8 Staphylococcal spp and the mecA gene. A blinded study of clinical blood culture samples compares Staph ID/R to an automated biochemical detection system.


Methods: Gene sequences were aligned and primers were identified that use Helicase-Dependent Amplification (HDA) to amplify conserved regions of the *tuf* and *mecA* genes. An array of hybridization probes was immobilized onto a silicon surface. For blood culture testing, a 50 uL blood culture aliquot was extracted enzymatically for 10 min, amplified for 30 min, and hybridized for 5 minutes producing produces eye-visible signals. Phenotypic determination of species results and methicillin-resistance was performed on the MicroScan. For discrepant results, definitive species identification was achieved by sequencing of *rpoB* and/or 16S rDNA genes, and definitive methicillin resistance was determined by growth on oxacillincontaining agar.


Results: The limit of detection for *S. aureus* was 3e3 CFU/mL, or 1 copy of input DNA, placing sensitivity on par with real-time PCR. 8 Staph spp, were identified and resolved with limits of detection ranging from 1-250 CFU/reaction. In a blinded study using clinical blood culture samples, the Microscan biochemical system correctly identified 93/98 samples that were mono-microbial infections for species and methicillin-resistance results, while Staph ID/R correctly identified 98/98 for bacterial species and methicillin-resistance. In 5 poly-microbial cultures Microscan was correct in 0/5 samples while Staph ID/R identified 4/5 samples correctly as defined by detection of all *Staphylococcal* species present in each sample.

Conclusion: Staph ID/R shows high sensitivity and specificity, equivalent to the gold standard for Staphyloccocal species identification, *rpoB* gene sequencing. Staph **19 Detect** ID/R represents an effective, low-cost, user-friendly approach to blood culture detection of *S aureus* and additional spp along with methicillin resistance. HDA and eyevisible detection together permit low-cost instrumentation. Clinical evaluations of the Portrait Dx platform are in progress at several sites.


Bacterial Identification Array










S. epidermidis









Capture probes are immobilized to the silicon surface. A diagonal set of control features verifies that chip orientation (COF), DNA extraction and amplification (EAC), and **detection** (HC, DC) functioned properly, validating test results.






Authors: B. Hicke¹, R. Jenison¹, C. Pasko¹, J. Dunn¹, H. Koss Jaeckel¹, D. Nieuwlandt¹, M. West¹, E. Woodruff¹, L. Rea¹, D. Weed², G. Denys³, X. Zheng⁴; 1. Great Basin Diagnostics, Salt Lake City, UT, 2. Denver, HIth. Med. Ctr., Denver, CO, 3. Indianapolis, IN, 4. Children's Mem. Hosp.. Northwestern Univ. Feinberg Sch. of Med., Chicago, IL.



CORPORATION



SA = NNNNNNNNNNANNANNASW = NNNNNNNNNNNTNNN

S. aureus

S. warneri

Limit of Detection

Bacterial		Limit of Detection (CFU/mL)		
Species	Gene	Pure	Mixed	
S. hominis/novo	tuf	1e3		
S. epidermidis	tuf	3e3	1e5	
S. capitis	tuf	1e5		
S. lugdunensis	tuf	1еб		
S. haemolyticus	tuf	1e4		
S. aureus	tuf	3e3	3e4	
S. warneri	tuf	1e5		
E. facaelis	tuf	ND		
MRSA/MRSE	mecA	1e3	1e3	

Blood culture samples were quantitated to determine CFU/mL, then serially diluted using a blood culture control. Staph ID/R was performed manually in 96-well format. *mecA* detection is more sensitive than species detection due to strength of capture probe for this sequence. LODs are on par with real-time PCR methods. LODs for mixed cultures were determined by placing a fixed amount of one competing bacterial species into a negative blood culture (3 x 10⁷ CFU/mL) and then titrating in the species of interest.

Great Basin Cartridge and Analyzer

The Staph ID/R assay is built into a injection-molded card. Reagents are lyophilized or placed in blister packs. In the instrument, optical sensors control motors that propel 10s to 100s of uL through channels and chambers. This mesofluidic-scale design and injection-molded plastic card, in combination with isothermal amplification and human eye-visible signal, enables a low-cost card and instrument.

The operator inserts ~50 uL blood culture into the sample port as shown, inserts the card into the desktop instrument, and initiates the test. Software automatically returns a result within 90 min. The report details the presence of staphylococci, specifically identifies the 12 Staph. species deemed most relevant, and indicates status of the drug resistance gene *mec*A.



Contact: rjenison@gbscience.com

Results

Limit of detection was determined for various Staphylococc strains using tuf (species) and mecA (methicillin-resistance) duplex amplification reactions, and ranged from 3 x 10³ to 1 x 10⁶ CFU/mL Analytical reactivity was demonstrated for 25 unique Staphylococcus strains representing 19 different species (data not shown). All amplified species reacted with the Staphylococcus genus control set except S. lentus and S. succinus, as expected. Each speciesspecific probe reacted with cognate *Staphylococcal* species amplicon with high specificity. The probe set described here showed no cross-reactivity with a panel of non-cognate bacteria except E. fa*caelis*, which reacted with the *Staphylococcus* genus probe set. This was expected and an *E. facaelis*-specific probe was created to distinguish it from *Staphylococcus* in the panel.

In a blinded study, 103 frozen blood culture samples positive for Gram positive cocci in clusters (GPCC) were tested. The distribution of species including methicillin-resistance determination is shown in Table 2. There was agreement on 93/98 samples from monomicrobial GPCC positive blood cultures between Staph ID/R and MicroScan. Resolution of discordant results with *RpoB* gene sequencing correlated in all cases with Staph ID/R.

Table 1.

Identification of Staphylococcus species: Micro Scan versus Portrait Staph ID/R

Resolved Result	Call by Test			
Species (n)	Portrait	Micro Scan		
S. aureus (19)	19	19		
MRSA (6)	6	6		
S. epidermidis (8)	8	8		
MRSE (40)	40	40		
S. hominis (2)	2	2		
MR S. hominis (4)	4	4		
MR S. haemolyticus (1)	1	0		
S. warneri/ pasteuri (1)	2	1		
MR S. warn/past (7)	6	6		
S. lugdunensis	0	0		
S. capitis/caprae (2)	2	1		
S. genus spp unknown (1)	1	0		
Non Staph E. facaelis (2)	2	2		
Non <i>Staph</i> (5)	5	3		
# Correct (%)	98 (100)	93 (94.9)		

For poly-microbial GPCC positive blood cultures, Staph ID/R found multiple Staphylococcus or E. facaelis organisms in each sample and was concordant with RpoB gene sequencing in all cases but one. MicroScan found only one Staphylococcus or E. facaelis organism in each sample. In two poly-microbial cases, MicroScan reported MRSA only whereas Staph ID/R detected a mixture of S. aureus, S. haemolyticus, and mecA in one and S. aureus, S. epidermidis, and mecA in another (Table 3). In one sample, mauve colonies were observed on ChromAgar MRSA, indicating the presence of MRSA. In the other case, mauve colonies were observed on ChromAgar SA, indicating S. aureus was present, but on ChromAgar only white colonies were observed. Sequencing of this colony indicated that the methcilillin-resistant organism was *S. epidermidis*.

Table 2. Analysis of Poly-microbial Specimens

	Test Results		Sub-culture Results				Sequencing Results
ID	MicroScan	Staph ID/R	TSA (unique colonies)	ChromAgar SA	ChromAgar MRSA	MRSA Screen	RpoB gene
10	MR S. epidermidis	S. epidermidis S. haemolyticus mecA	2	blue/green	blue/green	Mod	S. epidermidis S. haemolyticus
16	E. facaelis	E. facaelis S. epidermidis mecA	2	blue	blue	Mod	E. facaelis S. epidermidis
44	MRSA	S.aureus S. haemolyticus mecA	1	mauve	mauve and white	Light	S. haemolyticus S.aureus
58	MRSA	S. epidermidis S. aureus mecA	1	mauve	white	No	S. epidermidis S. aureus
102	E. facaelis K. pneumoniae	E. facaelis S. epidermidis S. hominis mecA	2	blue	blue and white	No	S. epidermidis S. cohnii E.facaelis

Conclusions

The Staph ID/R test was shown to rapidly identify the major Staphylococcus species and the mecA gene, which confers resistance to methicillin, directly from positive blood culture aliquots with high concordance compared to the standard of RpoB and 16S gene sequencing. Additionally, the test demonstrated the ability to detect more than one organism in a polymicrobial positive blood culture bottle.

Information about the species present and whether or not it harbors the mecA gene can aid in patient management. In the case of MRSA/MSSA tests the value has been clearly shown that a rapid test result lead to shorter hospital length of stay and lower treatment costs due to the initiation of appropriate therapy sooner.

By additionally providing knowledge about the potential methicillin-resistance status of coagulase-negative Staphylococcus, treatment can be adjusted to betalactam drugs for methicillin-sensitive infections which lead to improved outcomes compared with treatment with vancomycin, the generally accepted drug of choice.

Typically, ~30% of positive blood cultures are contaminants, not true infections. With the knowledge of species information in positive blood cultures, Staph ID/R could be used as an approach to distinguish true infections from contaminant blood cultures with high predictive value from just two positive culture bottles.